



3 Color

CD3 FITC-CD19 RPE-CD45 RPE-Cy5+

Product Exalpa's 3 color system, CD3 FITC-CD19 RPE-CD45 RPE-Cy5. Single laser excitation at 488 nm. All 3 color reagents are pre-titered and pre-tested for optimal performance.

Description Use Exalpa's 3 color CD markers to add flexibility to your current single laser system capabilities.

The CD3 epitope is expressed on the epsilon chain of the CD3/T cell antigen receptor (TcR) complex. CD3 is present on 65-85% of thymocytes and has a mitogenic effect on peripheral blood T cells. CD3 identifies human T cells expressing the 22-28,000 M.W surface antigen. CD19 identifies human B cells associated with 10% of peripheral blood lymphocytes. CD19 identifies lymphocytes expressing the 95,000 M.W. surface antigen. CD45 identifies all human leucocytes, lymphocytes, monocytes, polymorphonuclear cells, eosinphils and basophils in peripheral blood expressing the 180-220K M.W. T200 family surface antigen. CD45 is also present on thymus,spleen, tonsil and progenitor cells in bone marrow.

Isotype Mouse IgG1 kappa /IgG1 kappa/IgG1 kappa.

Clone M2AB (CD3 FITC)/1G9 (CD19 RPE)/ BHPT-6 (CD45 RPE-Cy5)

Fluorescence Information RPE-Cy5 is compatible with 488 nm single argon laser applications. With absorbence at 488nm and emission at 667 nm it is ideally suited for multi-color analysis on single laser instruments.

<u>Dye</u>	<u>Excitation</u>	<u>Emission</u>
FITC	488 nm	530 EM λ max nm
RPE	488 nm	565 EM λ max nm
RPE-Cy5	488 nm	667 EM λ max nm

Applications When using Exalpa's 3 color system CD3/19/45 it is appropriate to gate on CD45 RPE-Cy5 cell populations. Subsequent analysis may be done on the CD3 FITC and CD19 RPE cell populations.
Monitoring of T cells subsets in peripheral blood.
Characterization of subtypes of T cell leukemia's and lymphomas.
Studies of AIDS/HIV virus infection.
Analysis of CD3 complex related to the T cell antigen receptor.
Monitoring of B cells in peripheral blood.
Analysis of B cell subsets.
Study of B cell activation.
Study of B cell neoplasms.
Monitoring of Leukocyte cells subsets in peripheral blood.
Studies of progenitor cells in bone marrow.

For research use only. Not for use in human diagnostics or therapeutics.
Use PBMC: Add 10 µl of Exalph's 3 color antibody reagent/10⁶ PBMC in Exalpa Biologicals, Inc. Mix gently and incubate for 30 minutes at 2° to 8° C. Tel: 800.395.1137 or 978.425.1370 Fax: 866.924.5106 or 978.425.1376 Web: www.exalpa.com
Wash twice with PBS or fix with 0.5% v/v of paraformaldehyde in

PBS and analyze.

WARNING: Reagents containing sodium azide under acid conditions yields hydrazoic acid, an extremely toxic compound. Azide containing compounds should be discarded with proper caution. There are no warranties, expressed or implied, which extend beyond the description on the label of the product. Exalpha is not liable for property damage, personal injury or economic loss caused by the product. Nothing disclosed herein is to be construed as a recommendation to use Exalpha's products in violation of any patents. Exalpha cannot be responsible for patent infringement or other violations that may occur by the use of its products.

WHOLE BLOOD: Add 10 μ l of Exalpa's 3 color antibody reagent MAB/100 μ l of whole blood. Mix gently and incubate for 15 minutes at room temperature 20° C. Lyse the whole blood. Wash once with PBS or fix with 0.5% v/v of paraformaldehyde in PBS and analyze. See instrument manufacturers instructions for Lysed Whole Blood and Immunofluorescence analysis with a flow cytometer or microscope.

Storage Exalpa's 3 color reagents are supplied in PBS, 0.08% sodium azide and 0.2% carrier protein. Our conjugated monoclonal antibodies should be stored at 4 - 8° C until use. Exalpa's conjugated monoclonal antibodies should never be frozen. When stored properly they are stable for the period shown on the vial label.

Ordering Information	3 Color Form	Vial Size	Catalog#
	CD3/19/45	200 Test	X1066
	CD3/19/45	50 Test	X1066s

+ Portions of this product are manufactured under license from Carnegie Mellon University, U.S. Patent no. 5268486

For research use only. Not for use in human diagnostics or therapeutics.

For research use only. Not for use in human diagnostics or therapeutics.

Exalpa Biologicals, Inc., 2 Shaker Rd, #B101, Shirley, MA 01464
Tel: 800.395.1137 or 978.425.1370, Fax: 866.924.5100 or 978.425.1376, Web:
www.exalpha.com