



Bi-Test™ CD3 FITC - CD16 PE

Product: Anti-human T cell monoclonal I antibody and anti-human NK cell monoclonal antibody.

Description: The CD3 epitope is expressed on the epsilon chain of the CD3/T cell antigen receptor (TcR) complex. CD3 is expressed on 65-85% of thymocytes and has a mitogenic effect on peripheral blood T cells. Identification of human T cells expressing the 22-28,000 M. W. surface antigen. Identification of human NK cell antigen expressing the 50 - 70 kDa M.W. surface antigen associated with the IgG Fc receptor III on NK cells and Neutrophils. CD16 is expressed on approximately 15% of peripheral blood lymphocytes and is present on all resting NK cells. CD16 may be expressed on CD3 T cells from certain individuals.

Isotype: Mouse IgG1 kappa (FITC) and Mouse IgG1 kappa (PE)

Clones: M2AB (CD3 FITC) and J5511 (CD16 PE)

Applications: Monitoring of T cells subsets in peripheral blood; Characterization of subtypes of T cell leukemia's and lymphomas; Studies of AIDS/HIV virus infection; Analysis of CD3 complex related to the T cell antigen receptor; Monitoring of NK cell subsets in peripheral blood; Analysis of NK cell levels in peripheral blood; Characterization of subtypes of leukemias; Study of AIDS virus infection.

Use: PBMC: Add 10 μ l of MAB/10⁶ PBMC in 100 μ l PBS. Mix gently and incubate for 15 minutes at 2^o to 8^oC. Wash twice with PBS and analyze or fix with 0.5% v/v of paraformaldehyde in PBS and analyze.

WHOLE BLOOD: Add 10 μ l of MAB/100 μ l of whole blood. Mix gently and incubate for 15 minutes at room temperature 20^oC. Lyse the whole blood. Wash once with PBS and analyze or fix with 0.5% v/v of paraformaldehyde in PBS and analyze. See instrument manufacturer's instructions for Lysed Whole Blood and Immunofluorescence analysis with a flow cytometer or microscope.

Storage: Antibodies are supplied in PBS, 0.08% sodium azide and 0.2% protein carrier for FITC and PE. Antibodies should be stored at 4-8^oC. Monoclonal antibodies should not be frozen. Reagents are stable for the period shown on the vial label when stored properly.

Ordering Information:	Form	Vial Size	Catalog #
	Bi-Test™	100 Test	0316
	Bi-Test™	50 Test	0316S

For research use only. Not for use in human diagnostics or therapeutics.

REFERENCES:

1. Knowles RW. Immunochemical analysis of the T cell-specific antigens. In : Reinhert EL, Haynes BF, Nadl LM and Bernstein ID. eds. Leukocyte Typing II, Human T Lymphocytes. New York, NY: Springer-Verlag; 1986:259.
2. Kurrle R. Cluster Report:CD3. In:Knapp W, Dorken B, Gilks WR, Reiber EP, Schmidt RE, Stein H, and von dem Borne AEG Kr, eds. Leukocyte Typing IV, White cell Differentiation Antigens. Oxford, England: Oxford Press 1989:293.
3. Signal transduction via CD4,CD8 and CD28 in mature and immature thymocytes. Implications for thymic selection. Turka LA, Linsley PS, Paine R 3d, Schieven GI, Thompson GB, Ledbetter JA, J. Immunol. 1991 Mar :146(5): 1428-36.
4. T cell receptor/CD3-signaling induces death by apoptosis in human T cell receptor gamma delta + Tcells. Janssen O, Wesselborg S, Heckl-Ostreicher B, Pechhold K, Bender A, Schondelmaier S, Moldenhauer G, Kabelitz D I Immunol. 1991 Jan146(1):35-9.
5. Clonal analysis of human CD4-CD8-CD3- thymocytes highly purified from postnatal thymus Hori T, Spits H J. Immunol. 1991 Apr 146(7):2116-21.
6. Molecular cloning of the CD3 zeta subunit identifies a CD3 zeta-related product in thymus-derived cells Jin YJ, Claton LK, Howard FD, Koyasu S, Sieh M, Steinbrich R, Tarr GE, Reinherz EL. Proc Natl Acad Sci USA 1990 Ma: 87(9):3319-23.
7. Immunofluorescence Measurement in a Flow Cytometer using Low-Power Helium Neon Laser Excitation. Shapiro, H.M, Glazer, A.N., Christenson, L., Williams, J.M., and Strom, T. B. Cytometry 4,276, 1983.
8. Comparison of Helium Neon and Dye lasers for Excitation of Allophycocyanin. Loken, M.R., Kiej, J.F. and Kelly, K.,A. Cytometry 8, 96, 1987.
9. CD16 on human gamma delta T lymphocytes: expression, function, and specificity for mouse IgG isotypes. Braakman E; van de Winkel JG; van Krimpen BA; Jansze M; Bolhuis RL Cell Immunol 1992 Aug;143(1):97-107
10. Fc gamma RIII activation is different in CD16+ cytotoxic T lymphocytes and natural killer cells. Uciechowski P; Gessner JE; Schindler R; Schmidt RE Eur J Immunol 1992 Jun;22(6):1635-8
11. Involvement of a metalloprotease in spontaneous and phorbol ester-induced release of natural killer cell-associated Fc gamma RIII (CD16-II). Harrison D; Phillips JH; Lanier LL J Immunol 1991 Nov 15;147(10):3459-65
12. CD16- CD56+ natural killer cells after bone marrow transplantation. Jacobs R; Stoll M; Stratmann G; Leo R; Link H; Schmidt RE Blood 1992 Jun 15;79(12):3239-44
13. Natural killer function in flow cytometry: identification of human lymphoid subsets able to bind to the NK sensitive target K562. Vitale M; Zamai L; Neri LM; Manzoli L; Facchini A; Papa S Cytometry 1991;12(8):717-22
14. In vitro responsiveness to interleukins and theophylline of CD16+, CD56- natural killer cells in a patient with chronic granular lymphocyte disorder. Bayle C; Vitte-Mony I; Lang P; Pico J; Hercend T; Bertoglio J Leukemia 1992 May;6(5):470-6
15. Signal transduction by Fc gamma RIII (CD16) is mediated through the gamma chain. Wirthmueller U; Kurosaki T; Murakami MS; Ravetch JV J Exp Med 1992 May 1;175(5):1381-90

For research use only. Not for use in human diagnostics or therapeutics.

Exalpha Biologicals, Inc., 2 Shaker Rd, #B101, Shirley, MA 01464
Tel: 800.395.1137 or 978.425.1370, Fax: 866.924.5100 or 978.425.1376, Web: www.exalpha.com