



Exalpha Biologicals, Inc.

Human Placenta Tissue Control Lysate

CAT. No.
X1635C
SIZE

100 μ l

CONCENTRATION

1 mg/ml

APPLICATIONS

Western Blot

COMMENTS

Human placenta tissue lysate is ready to load on SDS-PAGE for Western blotting.

FORMULATION

Provided in 10% glycerol, 0.063 M Tris-HCl (pH 6.8), 2% SDS and 0.002% bromophenol blue, 5% 2-mercaptoethanol

STORAGE CUSTOMER

Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles

STABILITY

Products are stable for one year from purchase when stored properly

USE

Use 20 μ l of human placenta tissue lysate per lane. Lysate should be boiled prior to loading to ensure total denaturation of lysate. Boil for 10 minutes prior to loading.

For research use only. Not for use in human diagnostics or therapeutics.

Exalpha Biologicals, Inc. 5 Clock Tower Place, Suite 255 Maynard, MA 01754
Tel: 800.395.1137 Fax: 866.924.5100 www.exalpha.com info@exalpha.com

LAST MODIFIED
4/4/2007