



Human Placenta Tissue Control Lysate

BACKGROUND

ORDERING INFORMATION

CATALOG NUMBER

X1635C

SIZE

100 μ g

CUSTOMER STORAGE

Product should be stored at -20°C.
Aliquot to avoid freeze/thaw cycles

FORMULATION

Provided in 10% glycerol, 0.063 M Tris-HCl (pH 6.8), 2% SDS and 0.002% bromophenol blue, 5% 2-mercaptoethanol

SHIP CONDITIONS**STABILITY**

Products are stable for one year from purchase when stored properly

COMMENTS

Human placenta tissue lysate is ready to load on SDS-PAGE for Western blotting.

INSTRUCTIONS

Use 20 μ l of human placenta tissue lysate per lane. Lysate should be boiled prior to loading to ensure total denaturation of lysate. Boil for 10 minutes prior to loading.